# Codubix®

## GB INSTRUCTIONS FOR USE Prosthesis for cranioplasty

**PRODUCT CHARACTERISTICS**  
Codubix® cranial bone grafts are knitted of:


- polyester and polypropylene yarn – supplying sufficient strength against kinking and compression, as well as surface spread,
- polypropylene yarn – featuring low specific weight and low melting temperature allowing for appropriate stiffness and hardness of the graft.

High endurance, low specific weight, zero hydrophilic fluid absorbance, non-toxicity, chemical inactivity, high degree of healing in make Codubix® grafts easy to identify with the natural cranial bone due to their physical properties. Codubix® graft can be cut in order to shape it to fit the skull defect and it can be diagnosed without any artifacts.

**INDICATIONS**  
Codubix® knitted cranial bone grafts are intended for the reconstruction of skull cap defects.

**CONTRAINDICATIONS**  
No case of local intolerance to the implanted material has been reported so far. Codubix® knitted cranial bone grafts should be not used in children or where their use is limited by the process of physiological growth.

**GUIDELINES**  
It is critical to plan the skin incision line before the cranioplasty so that it goes at a minimum distance of 1 cm from the edge of the bone defect. When the skin flap has been peeled the periosteum should be cut along the edge of the bone defect and relocated as much outside as necessary to determine the edge of the bone defect. Both the rim of the defect and Codubix® graft are determined according to the given figures. The outer edge in Fig. 1 is bigger than the inner edge. The graft fits the defect as a wedge. Fig. 2 in order to fix Codubix® graft the openings must be made 5 mm from the edge of the graft and the edge of the bone. The openings are made with either a 2 mm dental drill or a trepan.

Haemostasis should be performed prior to the graft fixation. Following the introduction of fixation sutures through the openings resulting knots should be placed inside the openings. Correct adjustment of the graft as well as good point of contact of the graft edge with edge of the bone and tightening of the sutures safeguards the immobility of the graft in relation to the bone. The existing periosteum is tied to the graft with an atraumatic suture, thus covering the point of contact. In extensive cranioplasty it is essential to fix a sucking drainage in the area beneath the skin flap for 24 hours following the surgical procedure.

**APPLICATION**  
Prior to Codubix® graft implantation the following must be specified:

- size of the defect,
- shape of the defect,
- curvature relevant to that of the skull.

Final match and adjustment of the graft to the defect occurs intraoperatively (cutting of the graft is possible). After it has been modelled the graft needs to be carefully cleaned of the dust accumulated during the adjustment by rinsing in normal saline solution.

*\*) Author: Prof. PhD. Brzezinski, Dr. J. Gawlik, Neurosurgery Clinic, Neuropathy and Sense Organ Institute, Medical University of Łódź.*

The procedure of cleaning (time, volume of saline solution, etc.) may differ depending on the shape and size of pre-shaped prosthesis. It must be ensured that all visible particles are carefully cleaned.

**NOTE**

- Codubix® grafts are non sterile.
- Sterilisation should be performed using validated methods in an autoclave according to the parameters defined in the point sterilization.
- Minimum two-time sterilisation is acceptable if the prosthesis was not implanted before.
- Codubix® grafts must not be implanted after the expiry date.
- The product is intended for single use only. Due to the implantation method there is no possibility to reuse the product.
- Shaping, cleaning and re-packing of the graft should be performed under aseptic conditions.
- The product should not be steam sterilised in the manufacturer's packaging.

**STERILIZATION**  
Codubix® grafts are non sterile.  
Steam sterilisation in an autoclave is recommended using validated methods:

- temperature 121°C,
- 20 minutes.

It is recommended that the packaging protects the product against mechanical damage and complies with the requirements for steam sterilisation medical packaging.  
Personnel should have necessary competence for leading the steam sterilisation process (according to the Standard EN ISO 17665-1).

**COMPLICATIONS**  
Implantation of all cranial grafts may cause numerous kinds of complications, such as:

- in early post operative period: infection, haematoma, oedema of the skin flap, formation of haemorrhagic or serous exudate, increase in body temperature, ischaemia of the skin flap,
- in long post operative period: pains in the scar, sensitivity disorders, collapsing of the prosthesis.

All above mentioned complications may cause longer recovery, repeated surgical intervention or even re-implantation in extreme cases.

**STORAGE**  
Codubix® grafts should be stored:

- in dark rooms,
- at temperature range 15°C - 35°C,
- humidity limitation 25 - 65%,
- under conditions preventing from wetting and mechanical or chemical damage,
- in the original unit packaging (carton and direct packaging). Any mechanical damage to the packaging can expose the product to sterility loss.

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