

Dimensions of the final products in [cm]

- Female procedure gowns

Parameter / height	Sizes					Tolerances
	S	M	L	XL	XXL	
Chest circumference	100	108	116	128	134	± 2
Long sleeve length						± 1
158	57	57	57	57	57	
164	57	57	57	57	57	
170	60	60	60	60	60	
176	60	60	60	60	60	
Short sleeve length	24	24	24	24	24	± 1
Length of the gown						± 2
158	111	111	111	111	111	
164	111	111	111	111	111	
170	117	117	117	117	117	
176	117	117	117	117	117	