

































MATOSSET PICTOGRAMS - INSTRUCTIONS

How to read out the parameters given in one line:





Explanation of particular parameters:


Pictogram	Explanation	Pictogram	Explanation
Quantity		Product type	
1 (2 x 10)	1 piece tied: 2 piles by 10 pieces		drape: with adhesive tape, with opening, with opening and adhesive tape, under buttocks
	swab, swabs with Y and O cut neurosurgical swabs		T-shaped with an opening, for ophthalmology, for angiography
	tampons: ball, bean, whisker, triangle, thimble, grain, pepper		for urology, with U-shaped cut-out
	ribbon on a roll		for laparoscopy, for craniotomy, for cystoscopy
	abdominal swab with a cotton loop		vertical drape, Mayo table cover
	gynaecological towel Absorgyn		instrument table cover, 5-chamber pouch, Seni Soft absorbent underpad, film pouch for organs
	Absorba: absorbent swab, high-absorbent swab, tamponade dressing, cotton-wool and gauze tampon		wound closure strips, Medisorb dressing, tourniquet
	Absorba: eye dressing, nose swab		adhesive tape, on a bobbin pouch, medical fixing tape
	Tubula tubular bandage, bandage with a clip, Matosoft orthopaedic padding		
	adhesive tape, on a bobbin		

Pictogram	Explanation	Pictogram	Explanation
Product type		Product type	
	hand towel, nonwoven adhesive dressing Fixopore S		blade no. 10, no. 11, no. 15 blade no. 20, no. 24
	leg cover, arm cover, adhesive cannula fixation dressings		stitch cutter, scalpel
			anatomical forceps, umbilical cord clamp, surgical forceps
	medical equipment cover, C-arm cover, Camera Cover for medical device cables		syringe, needle, drain
	surgical gown: Standard, Perfect, apron sleeve protector, tied shoe cover, cap		hydrodelineation cannula 0.5x22mm, olive cannula 0.6x22mm, curved irrigation cannula
	surses cap, surgical cap with elastic band, with sweatband, with ties		Foley catheter, urine bag
	examination glove, surgical glove, face mask with ties, with elastic band		stitches, connector, 3-way stopcock with tubing
	container: square, round, tray, lamp cover		transfusion apparatus, venflon, non-return valve
	container: 2-chamber, 3-chamber, kidney dish		Redon bottle, monopolar handle, Codan line
	container: 3-chamber, for cables 2500ml, safe box for scalpels, safe box for scalpels and needles		Scalpfix scalp clip cartridge, 3-way stopcock
	sponge stick, cotton swab, tongue depressor		barraquer, drape clamp, neurosurgical swab
	eye spear for ophthalmology, plastic sharp-blunt scissors, Littauer scissors		Introcan safety needle
	scissors: sharp-blunt, sharp-sharp, for episiotomy		marker, baby measuring band, electrode cleaner
	scissors, umbilical scissors, clamps		3-way ramp of stopcocks
	sponge forceps, Kocher forceps/clamps, needle holder		
	Halsted Mosquito forceps, Pean forceps straight, Backhouse clamp umbilical clamp cur, skin staples remover		

Pictogram	Explanation
-----------	-------------

Symbols

13, 17, 20 
 30 g/m², 40 g/m² 



TMS, BP, Special, Comfort, TF,
 FB, PE



powder-free
 powdered

17 ≅ , 4 ≅

thread count for gauze
 basis weight of nonwoven
 gauze/nonwoven with X-ray detectable thread
 symbols of nonwoven/laminate/film types
 information on latex/powder presence/absence
 number of plies for gauze/nonwoven products

Pictogram	Explanation
-----------	-------------

Size

320x180 (Ø8)
 1 m x 2 cm
 200 ml
 nr 11
 M, L, XL, XXL
 6; 6,5; 7; 7,5; 8; 8,5; 9
 XS, S, M, L, XL

320 cm x 180 cm, opening of 8 cm diameter
 (the same units of measure)
 1 m x 2 cm (different units of measure)
 capacity 200 ml
 numbering used for describing blades
 sizes of surgical gowns
 sizes of surgical gloves
 sizes of examination gloves

Pictogram	Explanation
-----------	-------------