

**Patient gown - dimensions and tolerances given in [cm]**

<b>Parameter/size</b>	<b>Universal</b>	<b>height: 54/98</b>	<b>height: 104/128</b>	<b>height: 132/146</b>	<b>Tolerance</b>	<b>Test method</b>
Length of the gown from the highest shoulder point	110	50	70	90	± 2.0	EN 1773
Sleeve length	27	14	16	18	± 1.0	
Arm length	18	10	13	14	± 1.0	
Gown width at the bottom (unfolded)	140	92	112	120	± 2.0	
Length and width of the ties to tie at the waist	90 x 0.7				length ± 5.0 width ± 0.2	