

PRODUCT DIMENSIONS IN [CM]

Surgical clothing sets

Parameter/Height		Size						Tolerance	Test method	
		XS	S	M	L	XL	XXL			3XL
Chest circumference		96	104	110	116	120	124	130	± 2	EN 1773
Sleeve length		22						± 1		
Hip circumference		96	104	110	116	120	124	130	± 2	
Waist circumference (trousers)		90	95	100	105	110	115	120	± 2	
Top length	158-164	65	65	65	65	65	65	65	± 1	
	170-176	70	70	70	70	70	70	70		
	182-188	75	75	75	75	75	75	75		
Trousers length with belt	158-164	102	102	102	102	102	102	102	± 1	
	170-176	106	106	106	106	106	106	106		
	182-188	110	110	110	110	110	110	110		
Skirt length with belt	158-164	60	60	60	60	60	60	60	± 1	
	170-176	62	62	62	62	62	62	62		
	182-188	64	64	64	64	64	64	64		